

# USER'S MANUAL

## XREAL massage chair

# Dr.Care®



[www.drcare.us](http://www.drcare.us)



Made By  
**Dr.Care International Inc.**  
6565 E. Washington Blvd., Commerce,  
California, 90040, The U.S.A  
Tel: +1 (229) 800.0009  
[www.drcare.us](http://www.drcare.us)



# CONTENTS

- Introduction:..... 1
- Safety Introduction:.....2
- Product specifications:.....5
- Installation guide:.....6
- Features:.....15
- Name and function of each component:.....16
- Instructions of quick access buttons:.....17
- Pad controller guide:.....19
- Maintenance:.....45
- Trouble shooting:.....46
- Illustration:.....47
- Maintenance History:.....49
- Contents:.....52

# INTRODUCTION

*Thank you for choosing Dr.Care`s products!*

*XREAL massage chair DR-XR 929S is a high technology item to improve your health, to relax your body after a hard working day. There will be no more stress, any aches or pains, ... XREAL massage chair will bring you more for a happy life, an optimism, sound mind and sound body.*

*With over 24 years experience in the field of researching and manufacturing health care products, Dr.Care are non-stop researching and developing to create the modern unique leading products, more and more professional massage technique to satisfy Dr.Care`s customers.*

*To achieve the highest efficiency in treating pains and health recovery, please carefully read this User Manual book.*

*We always look forward to serving you better.*

*Your sincerely,*



Brian Williams

## SAFETY INTRODUCTION

- ▶ Precautions listed below will help to ensure proper use of this appliance.
- ▶ The hazard that may cause by incorrect use are classified into "WARNING" and "CAUTION", please ensure that you read all instructions before using the massage chair.

 <b>WARNING</b>	Improper operation may cause severe personal injury or substantial property damage.
 <b>CAUTION</b>	Improper operation may cause severe personal injury or productt damage.

	Prohibited		No Disassemble
--	------------	---	----------------

### Not Suitable For Below Persons

We recommend that those with any of the following conditions should consult your physician before using.

 <b>WARNING</b>	<ul style="list-style-type: none"> <li>◆ Those undergoing treatment, especially those who are not feeling well.</li> <li>◆ Those with malignant tumors or malignant abscesses.(Need to be supervised)</li> <li>◆ Women who are pregnant or menstruating.</li> <li>◆ Those with osteoporosis or fracture spine.</li> <li>◆ Those with skin disease or those whose skins are injured.</li> <li>◆ Those with high fever.</li> <li>◆ Close supervision is necessary when this product is used by, on, or near children, invalids, or disabled persons.</li> <li>◆ The appliance has a heated surface. Persins insensitive to heat must be careful when using the appliance.</li> <li>◆ This appliance can be used by children aged from 8 years and above and persons with reduced physical, sensory or mental capabilities or lack of experience and knowledge if they have been given supervision or instruction concerning use of the appliance in a safe way and understand the hazards involved.</li> <li>◆ Children shall not play with the appliance.</li> <li>◆ Cleaning and user maintenance shall not be made by children without supervision.</li> </ul>
--	---

### Controller Precautions

 <b>CAUTION</b>	<ul style="list-style-type: none"> <li>◆ Do not operate the controller or pull the controller cable with damp hand.</li> <li>◆ Do not spill water or other liquid into the controller to avoid fault.</li> <li>◆ Do not lay heavy objects on the controller.</li> </ul>
--	---

### Cautions Regarding the Adapters and Cord

## MAINTENANCE HISTORY

Date:	.....
Problems:	.....
Solution:	.....
Signed by Technician	
Date:	.....
Problems:	.....
Solution:	.....
Signed by Technician	
Date:	.....
Problems:	.....
Solution:	.....
Signed by Technician	

# MAINTENANCE HISTORY

Date:	.....
Problems:	.....
Solution:	..... <div style="text-align: right;">Signed by Technician</div>
Date:	.....
Problems:	.....
Solution:	..... <div style="text-align: right;">Signed by Technician</div>
Date:	.....
Problems:	.....
Solution:	..... <div style="text-align: right;">Signed by Technician</div>

# SAFETY INTRODUCTION

  
 CAUTION

- ◆ Always unplug this appliance from the electrical outlet immediately after using and before cleaning.
- ◆ An appliance should never be left unattended when plugged in. Unplug from outlet when not in use.
- ◆ Unplug from outlet during thinder or lightening storms, and before putting on or taking off parts.
- ◆ Never operate this appliance if it has a damaged cord or plug, if it is not working properly, if it has been dropped or damaged, or dropped into water the appliance to a service canter forexamination and repair.
- ◆ Check that the power coltage stated on the rating label corresponds to your electricity supply.
- ◆ Do not use a damp hand to unpulg the cord to avoid electric shock.
- ◆ Do not drag the cord, reach to the front plastic and to unplug it.
- ◆ Do not squeeze, over bend, pull hard, twist the cord.
- ◆ Keep the cord away from heated surfaces.
- ◆ Do not carry this appliance by supply cord or use cord as a handle.
- ◆ To discnnect, turn all controls to the off position, then remove plug from outlet.
- ◆ The appliance is only to be used with the power supply unit provided with the appliance.
- ◆ It must only be supplied at safety extra low voltage corresponding to the marking on the appliance.
- ◆ If the supply cord is damaged, it must be replaced by the manufacturer or its service agent or a simiarly qualified person in other to avoid a hazard.

## Best Circumstances

  
 CAUTION

- ◆ Do not p[er]erate under blanket or pillow. Excessive heating can occur and cause fire, electric shock, or injury to persons. Do not use this appliance where temperature is over 40 C.
- ◆ Do not use this appliance in humid or dusty environment, or when any part of the body is in contacts with plumbing or any similar ground.
- ◆ Never operate the appliance with the air openings blocked. Keep the air openings free of lint, hair and the like.
- ◆ Do not use outdoor.
- ◆ Please place the product on the flat level floor.
- ◆ Do not store near heat or open flame.
- ◆ Connect this unit to a properly grounded outlet only. See Grounding instructions.

# SAFETY INTRODUCTION

## Best Practices For Safe Use

 <b>CAUTION</b>	<ul style="list-style-type: none"> <li>◆ Only for users weighs below 120kg.</li> <li>◆ Use this product only for its intended use as described in this manual.</li>   <li>◆ Do not use the product for therapeutic treatment purpose.</li> <li>◆ Appliance should only be used by one person, do not use by two or more persons at the same time to avoid damage to the products.</li> <li>◆ Do not use the product right after meal, it may cause discomfort.</li> <li>◆ Remove the accessories from the cushions and other spots, if any damages are present, do not use it and refer to a qualified service professional for servicing.</li> <li>◆ If you feel any abnormal discomfort during the operation of this product immediately and do not use the product.</li> <li>◆ Use the appliance for 20 minutes each time is highly recommend</li> <li>◆ Make sure there is no obstacles behind the chair and there is enough clearance before adjusting the backrest.</li> <li>◆ Make sure there is no children, pet or other objects under the backrest or legrest. please be seated properly while using the appliance.</li> <li>◆ Do not drag or push the chair or roll on its caster on wood floor or carpet to avoid damage.</li> </ul>
 <b>PROHIBITED</b>	<ul style="list-style-type: none"> <li>◆ Do not stand on or stack heavy objects on the footrest.</li> <li>◆ Do not fall into sleep while the product is operating.</li> <li>◆ Avoid using the product under the influence of alcohol.</li> <li>◆ Never drop or insert any object into any openings.</li> <li>◆ If there is any smoke or unusual smell while using, stop using immediately, unplug the power cord.</li> </ul>

 <b>NO DISASSEMBLY</b>	<ul style="list-style-type: none"> <li>◆ Do not disassemble this product yourself to avoid failure.</li> <li>◆ Do not dismantle, fix or modify the yourself, please contact authorized dealer or service agent.</li> </ul>
---	--

# MAINTENANCE HISTORY




Date: .....	
Problems: .....	
Solution: .....	Signed by Technician
Date: .....	
Problems: .....	
Solution: .....	Signed by Technician
Date: .....	
Problems: .....	
Solution: .....	Signed by Technician

# ILLUSTRATION

# PRODUCT SPECIFICATIONS

## Zone Therapy Charts By Symptoms

Every Vertebra is related with the health, for example:

<p><b>Headache-Insomnia</b></p> 	<p><b>Shoulder</b></p> 	<p><b>Nerve Problem</b></p> 	<p><b>Knee Pain</b></p> 
<p><b>Heart Problem</b></p> 	<p><b>Diabetes</b></p> 	<p><b>Liver Problem</b></p> 	<p><b>Kidney Problem</b></p> 
<p><b>Rheumatism</b></p> 	<p><b>Arthritis</b></p> 	<p><b>Back Pain</b></p> 	<p><b>Sciatica</b></p> 

## XREAL massage chair

**PRODUCT NAME** : XREAL massage chair  
**MODEL** : DR-XR 929S  
**RATED VOLTAGE** : 120V~  
**RATED INPUT** : 230W 1.9A  
**RATED FREQUENCY** : 60Hz  
**SAFETY DESIGN** : CLASS I  
**SERIAL NUMBER** :  
 DR-XR 929S

**Dr.Care**  
better a little each day




Made by  
**Dr.Care International INC.**  
 6565 E. Washington Blvd., Commerce,  
 California, 90040, The U.S.A  
[www.drcare.us](http://www.drcare.us)

## INSTALLATION GUIDE

### Leg Extends/Retracts

Press and hold “” “” button to extend or retract the footrest to the length you desire for better massage experience, release the button to fix the position.



A footrest sensor is set at the bottom of footrest, the footrest will adjust the length automatically once the chair started.

### WARNING

- ▶ Do not play on the massage chair or operate it improperly.
- ▶ Ensure there is no obstruction when raising or reclining the backrest and footrest.
- ▶ Do not place foot, fingers or any other objects in the space between the seat and the backrest/footrest when it is in the reclining or raising motion.
- ▶ Ensure massage chair is restored back to its original position before getting off the massage chair.


### After Use

Always unplug the massage chair from the electrical outlet immediately after use. Secure the Touchscreen Controller onto the holder. Turn off the power switch and unplug the massage chair from the electrical outlet while not in use to avoid children power on the massage chair accidentally. When not in use for a long period of time, we recommend that the chair be covered, the power cord be coiled and that the unit be placed in an environment free of dust and moisture.

## ILLUSTRATION

### Acupressure Point Massage Instruction

Every Vertebra is related with the health, for example:



- |   |   |
|---|---|
| 1. The first bone of neck vertebra diagnoses headache.                      | 6. The ninth bone of chest vertebra diagnoses chololith.                  |
| 2. The third bone of neck vertebra diagnoses hard of hearing.               | 7. The sixth bone of chest vertebra diagnoses nephropathy.                |
| 3. The fifth bone of neck vertebra diagnoses bronchus.                      | 8. The third bone of waist vertebra diagnoses constipation and back pain. |
| 4. The second bone of chest vertebra diagnoses cardiopathy.                 | 9. The fifth bone of waist vertebra diagnoses hemorrhoid.                 |
| 5. The fourth bone of chest vertebra diagnoses liver disease and tummy bug. | 10. Hemorrhoids, anus, enhance physiology.                                |

## TROUBLE SHOOTING

No.	Symptoms	Possible Causes	Possible Solutions
1	Noise is heard while in use	This sound is normal and comes from the motor or massager rubbing against the cloth cover	No need to take any measures
2	Appliance does not start	The power cord is not secured into the socket	Insert the power plug properly into the power socket
		Power switch is off	Turn on the switch
		The fuse is blown	Replace only with a fuse of same type and rating
3	Airbag malfunction	The air hose is blocked	Straight the air hose and clear out the blockage

If problem remains, please contact the dealer for maintenance.


Important: Do not attempt to repair this product yourself, our company shall not be responsible for the injury or damage caused by any servicing without authorization.

Thank you for your cooperation.

## INSTALLATION GUIDE

### Commissioning

- ▶ Use the controller to recline the backrest to its lowest position, raise and extend the footrest to its top position, make sure there is enough clearance before commissioning.
- ▶ Return the backrest and footrest back to its default position after commissioning and turn off the power switch.



Before Sitting in the Chair

- ▶ Ensure the footrest is well assembled before sitting on the massage chair.
- ▶ Ensure the massage chair is returned to its original position before sitting on it. If there is any part not in default position, please turn off the massage chair to allow it to restore.

#### CAUTION

- ▶ Do not take seat while the footrest is lifted.
- ▶ Do not put excessive load on the footrest to avoid damage.

#### DANGER

- ▶ Improper connection of the equipment grounding conductor can result in a risk of electric shock. Check with a qualified electrician or serviceman if you are in doubt as to whether the product is properly grounded. Do not modify the plug provided with the product – if it will not fit the outlet, have a proper outlet installed by a qualified electrician.


#### EARTHING

- ▶ This massage chair belongs to Level One electrical appliance, make sure to use three-core power supply socket and to guarantee the ground wire connected to the power supply socket has obtained good ground connection to avoid electric leakage, electric shock and some other negative effect during use.

## INSTALLATION GUIDE

### Moving the Massage Chair

- ▶ Roll the chair on its casters to the desired location.
- ▶ Incline the backrest to upright position.
- ▶ Turn off the power switch and unplug the appliance from the electrical outlet.
- ▶ Using two people, push down on the backrest and lift up on the front side by grabbing the moving Holder located to the front bottom of the armrest, roll the chair on its casters to the desired location (check the sketch for reference).




#### **⚠ CAUTION**

- ▶ Do not move the chair when there is user inside.
- ▶ To carry this product, do not lift by the footrest or armrest, hold the chassis and steady the backrest (have 4 or more people to carry the product).

### Connecting to Power

- ▶ Connect this appliance to a 3-wire properly grounded outlet.
- ▶ Turn on the main power switch by the right hand side of the massage chair.



#### **⚠ CAUTION**

- ▶ Please verify that there is no damage to the supply cord before switching on the massage chair.
- ▶ Please verify that the main power switch is turned off before connecting to power.

## MAINTENANCE

### Care and Maintenance Instructions

- ▶ To treat stains on back panel, armrest, footrest, plastic parts or air hose, spot clean with mild detergent foam, then wipe it with a dry cleaning cloth.
- ▶ To treat stains on controller, power box, spot clean with a dry cleaning cloth.
- ▶ To treat stains on back cushion or pillows, use a water-based upholstery cleaner only, wipe it with moist cloth, then air dry.
- ▶ To clean the foot/calf unit covers, remove it from the chair, hand wash with a mild detergent, then line dry.

#### **⚠ CAUTION**

- ▶ Pull out the power cord from the power strips before care and maintenance. To avoid electrical shock or injury, do not touch the power cord with wet hands. Make sure no detergent enter into the power box while cleaning.
- ▶ Do not use benzene, thinner or other solvents on your product, it may cause fading or damage to the product.
- ▶ Do not iron the upholstery.

### Storage Methods

- ▶ Keep the product free from dust. Cover the product with a dust-free cloth to prevent collection of dust if not in use of a long period of time.

#### **⚠ CAUTION**

- ▶ Do not leave this product exposed to direct sunlight for long periods of time. This may cause fading or damage.

### Warm Tips

During daily use, please check for the following symptoms:

- ▶ Burnt smell.
- ▶ The product experiences unusual power ON/OFF when touching the power cord.
- ▶ Hot power cord.
- ▶ Other abnormal symptoms.

#### **⚠ CAUTION**

- ▶ Please stop using the product to avoid breakdowns or accidents if any above symptom spotted.
- ▶ Switch off the power button and disconnect the plug, contact the local distributors or supplier for repair.
- ▶ All other servicing except above care and maintenance should be serviced by an authorized dealer or service center.

#### **⚠ WARNING**

Pull out the power cord from the power strips before care and maintenance. To avoid electrical shock or injury, do not handle the power cord with wet hands. Do not try to dismantle or repair the product yourself, please send your massage chair to the authorized service center.

## PAD CONTROLLER GUIDE

About Pad	
<b>Serial Number of Pad</b>	3c6617935691d978
<b>MAC Address</b>	28:ed:e0:65:58:2a
<b>Software Version</b>	1.0.42
<b>System Version</b>	v1.1.0
<b>Hardware Version</b>	v1.0
<b>User Authorization</b>	<b>Dr.Care</b>


Help	
<b>If you are experiencing difficulties, consult the user manual that came with your chair.</b>	
<b>Bluetooth Connection Troubleshooting</b>	
1. In your device settings, check to make sure that Bluetooth is switched. The Bluetooth icon shows on the home display screen if when it is on. Go to Settings > Bluetooth, and switch on Bluetooth to pair massage chair.	
2. Check to make sure the chair is connected to power.	
<b>If problem persists, please contact us.</b>	
<b>Business hours: Monday to Friday, 9AM - 5PM PST</b>	

## INSTALLATION GUIDE


### Position & Placement for Installation

Clearance space for the massage chair:

- Ensure there is an adequate clearance for the installation: It needs 10cm clearance for the reclination of the backrest; It needs 80cm clearance space for the extension of footrest and sliding forward of the chair.
- Keep it 1 meter away from the TV, radio or other audio&video appliance to avoid signal interference.



- Keep the distance between the power socket and power switch within 1.5 meter.




### ⚠ CAUTION

- Keep the distance between the power socket and power switch within 1.5 meter to avoid damage of power cord while massage chair sliding forward.
- Straighten the power cord while in use to avoid tangle or damage.
- When the massage chair sliding forwards or back, please check to make sure there are no children, pets or other obstacles around to avoid accident.

## INSTALLATION GUIDE

### Floor Protection

- We recommend to place a piece of carpet or soft mat underneath it to avoid damaging the floor.



### Best Circumstances



Do not use the massage chair in high moisture environment such as near swimming pool or bathroom to prevent electric shock accident.



Please keep the massage chair away from heated surfaces to prevent fire hazard or damage to the product.



### Grounding Instructions

- This product must be grounded. If it should malfunction or breakdown, grounding provides a path of least resistance for electric current to reduce the risk of electric shock. This product is equipped with a cord having an equipment - grounding conductor and a grounding plug. The plug must be plugged into an appropriate outlet that is properly installed and grounded in accordance with all local codes and ordinances.



This product must be grounded, the plug must be plugged into an appropriate 3-wire grounded outlet that is properly installed and grounded in accordance with all local codes and ordinances to reduce risk of electrical shock.

## PAD CONTROLLER GUIDE



## PAD CONTROLLER GUIDE

**CONNECT Massage Chair:** you can connect this pad controller with the massage chair to control all functions of the chair.

Choose model number of this massage chair DR-XR 929S, and connect it.

**UPGRADE massage chair:** This stage is specially used for engineer jobs, to load the program into the chair and to control the mechanism.

Please DO NOT connect or load this chair with any other programs (unless it was provided by Dr.Care's engineer).

**SOUNDS:** You can adjust volume of this pad controller as you wish.

**SCREEN BRIGHTNESS:** You can adjust the light brightness of this pad controller as you wish.

You also can adjust the screen timeout: 1 minute, 2 minutes and 3 minutes as you wish.

**HELP:** You can troubleshoot the Bluetooth Connection.

Please make sure that the chair is connected to power supply.

**USER GUIDE:** You can learn how to use / control the chair.


**LANGUAGE:** You can choose your language to control the chair.

In future we will have full list of many languages for your choice. Currently we have only English.

**ABOUT:** In this menu, you can review full information about this pad controller like: Serial number, software version, system version, hardware version, user authorization

**In this User Authorization:** There is Charging Permission tab.

You can turn it ON/OFF this charging function. When it is OFF, wireless charge and center console charging will stop working.



## INSTALLATION GUIDE


### Installation Method

#### Step 1, Take all massage chair units out of the cartons

- ▶ Unpack the cartons: take out all upholsteries and accessories out of the cartons.
- ▶ Before beginning installation, please verify that all components are included (check with the packing list carefully)

#### Step 2, Install the shoulder massage airbag and speaker unit


- ▶ Connect the shoulder airbag and speaker unit to the backrest connecting port.
- ▶ Line up the two connector pillar slots of shoulder airbag and speaker unit with the pillars on the massage backrest, press to insert into place, then secure with 3 screws(Fig. 2)  
(Same installation method with the other side)




## INSTALLATION GUIDE

### Step 3, Install Armrest Unit


- ▶ Connect the air valve, LED light connectors and air hoses of the armrest unit to the matching connection port and nozzles on the armrest frame( Follow the sequence shows in below figure);



- ▶ Line up the armrest unit groove with the armrest frame and press downwards, fit the screw to secure the armrest in place(a). (Same installation method with the other side)
- ▶ Put on the cap cover(tell the right from the left), then rotate clockwise to tighten it up(b).
- ▶ Fit a screw from inner side of armrest to secure(c).
- ▶ Insert the lower end of touchscreen controller holder into the socket, turn the nut clockwise to secure the holder into place.





## PAD CONTROLLER GUIDE



**CONNECTION:** It allows you to connect the pad controller with massage chair and upgrade the program of massage chair.

## PAD CONTROLLER GUIDE




You can check **Connection**, **Sounds**, **Screen Brightness**, **Help**, **Language**, **User guide** and **About** in Settings.


## INSTALLATION GUIDE

### Step 4, Install the Footrest Unit


- ▶ Thread the plastic bushing onto the footrest massage unit connection pin.



- ▶ Remove the screw off the coupler bracket, latch upwards to unlock it, fit the footrest unit connector to the coupler socket and latch downwards to lock it, tighten the screw.



- ▶ Raise the legrest unit, thread the air hose into the air nozzle, and fit legrest drive wire connector into the drive wire port, make sure they are well secured. (Check below illustration for reference)
- ▶ Zip up the cover and cushion.





## INSTALLATION GUIDE

## PAD CONTROLLER GUIDE

### Step 5, Other Accessories


Connect the power cable to the product power socket located at the bottom of right armrest.

Connect power to chair, switch ON the power button, use the controller to raise the backrest to fully upright position.



Settings

## PAD CONTROLLER GUIDE



**ADJUSTMENTS:** You can control the time of massage session and adjust the position of the chair, backrest and footrest.

**ZERO Gravity:** There are 2 positions of Zero Gravity: # 1 and # 2. Enjoy any position you like. You also can turn ON/OFF this Zero Gravity.

**MASSAGE TIME:** There are 3 options: 10 minutes, 20 minutes and 30 minutes. Please choose it.

**RECLINE:** Press the arrow UP and hold it to move the backrest up. Press the arrow DOWN and hold it to move the backrest down.


**LEGS:** Press the arrow UP and hold it to move the legs up. Press the arrow DOWN and hold it to move the legs down.

**FEET:** Press the arrow UP and hold it to move the feet massage parts up. Press the arrow DOWN and hold it to move the feet massage parts down.

## FEATURES

### 1. One Press Control Sliding

It will slide forward by 30cm and recline afterwards. Glider recliner linkage mechanism to save space and provide utmost relaxing position.



### 2. Zero Gravity Position

At such a weightless position that the height of the legs slightly higher than the upper body, the weight of body can evenly distributed throughout, reducing the pressure of the spine and joints both physically and mentally.

### 3. 4D Deep Fascia Massage

Innovative 4D double action deep fascia massage mechanism and L- shaped guiding track conform to ergonomic principle.

### 4. Auto/Manual Footrest Extension

The footrest can extend automatically to customize the user's leg length, you may also use the controller button to adjust the footrest according to your needs. The footrest can extend by 17cm automatically for various users of different height ranges.

### 5. Whole Body Air Pressure Massage

Whole body air pressure massage for shoulders, arms, back, waist, hips, legs and feet.

### 6. Rolling Massage

3 underfoot rollers provide all - direction pressure massage for the reflex zone.

### 7. Versatile Massage Techniques

Mimic human hand kneading, knocking, sync, tapping, shiatsu massage techniques with extra calf rubbing massage for custom - tailored relaxation.


### 8. Quick Control

Quick buttons at armrest meets the demands for easier, finger-tip control.


### 9. Super Long L-Style Track

## NAME AND FUNCTION OF EACH COMPONENT

### External Structure



### Internal Structure



## PAD CONTROLLER GUIDE

Press the Power/ sleep button to turn on the Touchscreen Controller (verify that it is connected to the massage chair).

**Swipe left** on Home screen to go to **Menu** screen.

Tap the **Adjustments** icon on Menu screen to go Adjustments screen.

You can adjust **Zero gravity, Massage Time, Recline, Legs and Feet.**


**Zero Gravity** Zero 1, Zero 2 and Off (turn off Zero Gravity).  
Tap the icon to start.

**Massage Time** 10mins, 20mins and 30mins



**Recline** Tap and hold the icon to recline/ incline the massage chair back.

**Legs** Tap and hold the icon to lift/ lower the massage chair footrest.

**Feet** Tap and hold the icon to extend/ retract the footrest.



## PAD CONTROLLER GUIDE



You can adjust **Zero gravity**, **Massage Time**, **Recline**, **Legs**, **Feet** in Adjustments.

## INSTRUCTIONS OF QUICK ACCESS BUTTONS

**Power/Auto Massage Button**

► Long press rotary switch to power on/off massage chair, quick press rotary switch to toggle among 6 auto massage programs.

**Rotary Switch**

► Rotate clockwise/counterclockwise to increase/decrease 3D massage intensity.

**Zero Gravity Position**

► Press to put the chair in Zero Gravity Position 1, press again to put the chair in Zero Gravity Position 2, press for the third time to restore to default position.

**Feet Extends**

► Press and hold to extend the footrest; release the button to stop. (Make sure there is no obstruction in front)

**Feet Retracts**

► Press and hold to retract the footrest; release the button to stop.

**Pause**

► Press to pause massage chair function. Press this button again to continue.

**Feet Up**

► Press and hold to raise the footrest; release the button to stop.

**Feet Down**

► Press and hold to lower the footrest, release the button to stop.

**Back Heat Button**

► Press to turn on back heat, press again to turn off back heat.

**Back Up**

► Incline: Press and hold to incline the backrest to upright position; release the button to stop.

**Back Down**

► Recline: Press and hold to recline the backrest; release the button to stop.



**Touchscreen Controller User Guide**

► See the "Touchscreen Controller User Guide" at touch screen controller display.

# INSTRUCTIONS OF QUICK ACCESS BUTTONS


# PAD CONTROLLER GUIDE

## Controller



Adjustments

## PAD CONTROLLER GUIDE



**AIRBAGS:** Turn ON/OFF airbags at these positions on your body: Full of body, shoulder, Back, Legs, Hips. Press once to turn ON, press once again to turn OFF.


**AIR STRENGTH:** You can choose and control the strength of air pressure on you body, there are 5 levels: weak, medium, strong, stronger and strongest.

**FOOT ROLLER:** you can turn ON/OFF the roller massage at your foot. You also can control the speed of massage, there are 3 levels of massage speed: weak, medium and strong.

**BACK HEAT:** You can control the temperature of heating. there are 3 levels of heating: low, medium and high. You also can turn ON/OFF the heating if you like.


**OPTIONS:** You can control airbags, foot roller, calf rub and heating.

## PAD CONTROLLER GUIDE




Home Screen

Home Screen may vary depending on the task you're performing.  
(Massage chair connection/ Recommendations/ Running massage program controls)  
Screen frame is highlighted when task is active.



No Massage Chair Connection  
Tap **Pair Chair** to connect massage chair.

## PAD CONTROLLER GUIDE



**DEMO:** Choose and hit DEMO to enjoy all special techniques of massage, enjoy the unique massage functions.

**AUTO:** In this Auto, there are 12 modes of automatic massage combinations, you can choose any Auto massage you like.

#### 12 AUTO MODES:

1.Stretch, 2.Reborn, 3.Tender Care, 4.Royal Treatment.  
5.Stretch Relief, 6.Gentle Care, 7.Neck & Shoulder, 8.Back & Spine.  
9.Waist & Hip, 10.Morning Buzz, 11.Refreshing Nap, 12.Sweet Dreams.

**MANUAL:** You can choose any technique of massage you like.

**OPTIONS:** You will see details of each function, choose and control it by yourself.

**ADJUSTMENTS:** You can adjust weak, medium, stronger, zero Gravity....you like.

**SETTINGS:** You can set up all about this pad controller.

## PAD CONTROLLER GUIDE

Press the Power/ sleep button to turn on the Touchscreen Controller (verify that it is connected to the massage chair).

Swipe left on Home screen to go to Menu screen.

Tap the Options icon on Menu screen to go Options screen.

You can choose **Airbags, Air Strength, Foot Rollers, Calf Rub and Back Heat**.

**Airbags** Kneading, Tapping, Shiatsu, Sync, Knocking and 3D Shiatsu.  
Tap the icon to start.

**Air Strength** 5 levels and Off ( turn off Air Strength).


**Foot Rollers** 3 levels and Off ( turn off foot Rollers).

This option is not adjustable if full or Legs Airbags is selected.



**Calf Rub** 3 levels and off ( turn off Calf Rub).

This option is not adjustable if Full or Legs Airbags is selected.

**Back Heat** 3 levels and off ( turn off 4D Strength).




# PAD CONTROLLER GUIDE



This page shows **Airbags, Air Strength, Foot Rollers, Calf Rub, Back Heat options**, you can adjust the Options setting during the massage session.

# PAD CONTROLLER GUIDE



Auto Massage Programs

## PAD CONTROLLER GUIDE

### Method 1


Press the Power/ sleep button to turn on the Touchscreen Controller (verify that it is connected to the massage chair).



There are **3 featured auto message programs** on Home screen.



**Tap** the program you like to start.



### Method 2

Press the Power/ sleep button to turn on the Touchscreen Controller (verify that it is connected to the massage chair).

Swipe left on Home screen to go Menu screen.

Tap the Auto icon on **Menu**  screen to go Auto screen.

You may choose among 12 pre-set Auto wellness massage program. Tap the “◀/▶” icon at both sides of the screen to scroll to other Auto massage programs.

Tap the program you like to start.

## PAD CONTROLLER GUIDE



Options

## PAD CONTROLLER GUIDE

**MANUAL:** ——— > you can play any functions of massage you like, control the chair, order the chair to serve you.

In this Manual, there are 6 techniques of massage: Kneading, Sync, Knocking, Tapping, Shiatsu, 3D Shiatsu.

**ADJUST:** This will help you to move up and down the back rollers when you choose massage at “a point” on your body.


**PART:** you can choose and activate the massage on any part of body. Whole body, Partial body and a point of body.

**WIDTH:** you can control the massage hands to move wide, wider, medium and narrow.



**3D STRENGTH:** There are 5 levels of 3D press deeper strength, you can control to increase or decrease the strength at any level you like. This 3D roller can press deeper on any part of your body to bring you the feeling of real hand from massage experts.

**4D STRENGTH:** There are 5 levels of 4D press deeper strength in the opposite dimension. You can control to increase or decrease the strength at any level you like. This 4D roller can press deeper on any part of your body to bring you the feeling of real hand from massage experts.


**SPEED:** You can choose and control the speed of robot roller (massage hands): Slow, medium, fast, faster and fastest. There are 5 levels for your choice.




## PAD CONTROLLER GUIDE



## PAD CONTROLLER GUIDE



### Reset



## PAD CONTROLLER GUIDE


Press the Power/ sleep button to turn on the Touchscreen Controller (verify that it is connected to the massage chair).

Swipe left on Home screen to go to Menu screen.



Tap the Manual icon on Menu screen to go Manual screen.

You can choose the massage Technique you like, and set the massage Part, Width, 3D Strength, 4D Strength and Speed.

- Techniques** Kneading, Tapping, Shiatsu, Sync, Knocking and 3D Shiatsu. Tap the icon to start.
- Part** Whole, Partial and point. When Point massage is working, Poin focus icon will show on the screen. Tap and hold the icon to adjust the massage focus you desire.
- Width** Wide, medium and Narrow.
- 3D Strength** 5 levels and off ( turn off 3D Strength).
- 4D Strength** 5 levels and off ( turn off 4D Strength).
- Speed** 6 levels. Speed is not adjustable if Shiatsu or 3D Shiatsu massage is selected. (⊘)




## PAD CONTROLLER GUIDE




Tap the massage technique and the way you like to start, you can adjust the **manual massage settings** during the massage session.


## PAD CONTROLLER GUIDE



Body scan




# PAD CONTROLLER GUIDE



## Massage Program Running

Home screen shows status of current massage session.  
Tap to stop/ pause the massage session.



Turn on this icon, the chair will reset to original position (stand up right) after finishing a massage session.

# PAD CONTROLLER GUIDE



## Manual Massage